

Electrodrive patient mover

COMPATIBILITY CHART

Electrodrive's powered bed movers are designed to be used with hospital beds, stretchers and trolleys.

This compatibility chart lists products that have been tested successfully with one or more of our patient movers and approved for use by Electrodrive.



Gzunda GZ30



Gzunda GZBARI (bariatric)



Gzunda GZS



Powered Fifth Wheel (P5W)

If you can't find the product you are looking for, please contact our Healthcare manager on 03 9300 8500.



COMPATIBILITY CHECK

Manufacturer	Model/name	GZ30	GZBARI	GZS	P5W	Accessory bracket	Notes
Acare	HCB 9433-H Trendelenburg and Reverse Electric Bed	✓	✓	✓	X		
Arjo	500	✓	✓	✓	✓		
Arjo	Contoura 380	✓	✓	✓	✓		
Arjo	Contoura 460	✓	✓	✓	✓		
Arjo	Contoura 480	✓	✓	✓	✓		
Arjo	Contoura 1080	✓	✓	✓	✓		Foot extension*
Arjo	Enterprise 5000	✓	✓	✓	✓		Foot extension*
Arjo	Enterprise 8000	✓	✓	✓	✓		Foot extension*
Arjo	Enterprise 9000 (without Indigo Powered 5th Wheel)	✓	✓	✓	X		Foot extension*
Arjo	Sorrento	✓	✓	✓	✓		
CGM Taiwan	Vanguard Classic	✓	✓	✓	✓	Recommend use with SIE0369 (GZ30SG 'Singapore Slimline')	
Engrit	EN71	✓	✓	✓	✓		Hitch from head-end only
Hausted	Unicare 800	X	X	✓	✓		
Horizon	Electric Powered Stretcher	X	X	✓	X		Column bed
Hausted	Youth Stretcher	X	X	✓	X		Column bed
Hendicare	700	✓	✓	✓	✓		
Hill-Rom	305 (Manual bed)	✓	✓	✓	✓		
Hill-Rom	405 (Electric bed)	✓	✓	✓	✓		
Hill-Rom	Avantgarde 800	X	X	✓	X		Column bed
Hill-Rom	Avantgarde 1200	X	X	✓	X		Column bed
Hill-Rom	Avantgarde 1600	X	X	✓	X		Column bed
Hill-Rom	Bariatric Electric Bed: Compella dynamic, Affinity 400	✓	✓	✓	✓		
Hill-Rom	Birthing Electric Bed: Affinity 600	✓	✓	✓	✓		
Hill-Rom	Centrella 355 Smart Bed	✓	✓	✓	✓		
Hill-Rom	Centuris	✓	✓	✓	X		
Hill-Rom	Emergency Electric Bed: HR900SR	✓	X	✓	✓		
Hill-Rom	FE5000	✓	✓	✓	✓		
Hill-Rom	HLE 8000	✓	✓	✓	✓		
Hill-Rom	HR900	✓	✓	✓	✓		Suitable with gripper
Hill-Rom	Medicraft 200 (MC200)	X	X	✓	X		
Hill-Rom	Medicraft 500 (MC500)	✓	✓	✓	✓		




COMPATIBILITY CHECK

Manufacturer	Model/name	GZ30	GZBARI	GZS	P5W	Accessory bracket	Notes
Hill-Rom	Medicraft 700 (MC700)	✓	✓	✓	✓		
Hill-Rom	Medicraft FES400 (Bariatric)	X	X	✓	X		
Hill-Rom	Standard Electric Bed: HR900 Accella	✓	X	✓	✓		
Hill-Rom	Total Care Connect	X	X	✓	X		Column bed
Hill-Rom	Transport Stretcher P8005	X	X	✓	X		
Hill-Rom	Versacare 3200	✓	✓	✓	✓		
Howard Wright	M7	✓	✓	✓	✓	If using gripper, must use Howard Wright gripper SIE0305 (jaw-lock, 50 mm)	Foot extension*
Howard Wright	M8	X	X	✓	X		Column bed
Howard Wright	M9	✓	✓	✓	✓		Foot extension*/requires 90 mm tine
Ikan	Rose Floorline	X	X	✓	X		Column bed
Invacare	1600	✓	✓	✓	✓		
Invacare	Alegio	✓	✓	✓	✓		
Invacare	Aria	X	X	✓	X		Column bed
Invacare	Etude Basic	✓	✓	✓	✓		
Invacare	Etude Classic	✓	✓	✓	✓		
Invacare	Etude Duo	✓	✓	✓	✓		
Invacare	Octave	✓	✓	✓	✓		
Invacare	ScanBed 750	X	X	✓	X		
Liftcare	Protean	X	X	✓	X		Column bed
Linet	Eleganza Smart	X	X	✓	X		Column bed
Linet	Eleganza Standard	X	X	✓	X		Column bed
Linet	Eleganza XC	X	X	✓	X		Column bed
Jackson Shand	ECO	✓	✓	✓	✓		
Jackson Shand	Phoenix	✓	✓	✓	✓		
Jomor	FLO101	X	X	✓	X		Column bed
Jomor	LHL272	✓	✓	✓	✓		
Jomor	ULT222	✓	✓	✓	✓		
Joyce	400/E	✓	✓	✓	✓		
Joyce	900/E	✓	✓	✓	✓		
Joyce	Centenary	✓	✓	✓	✓		
Maxi Promotions	Deluxe Regal	✓	✓	✓	✓		
Maxi Promotions	Floor Care	X	X	✓	X		Column bed
Maxi Promotions	Low Level	✓	✓	✓	✓		
Maxi Promotions	Ultimate Regal	✓	✓	✓	✓		



COMPATIBILITY CHECK

Manufacturer	Model/name	GZ30	GZBARI	GZS	P5W	Accessory bracket	Notes
Modsel	Contour Classic	✓	✓	✓	X		Use SIE0412 gripper bar with extension
Modsel	Contour Endosurge	✓	✓	✓	X		Use SIE0412 gripper bar with extension
Modsel	Contour Multi-X	✓	✓	✓	X		Use SIE0412 gripper bar with extension
Modsel	Contour Portare	✓	✓	✓	X		Use SIE0412 gripper bar with extension
Modsel	Contour Portare-X	✓	✓	✓	X		Use SIE0412 gripper bar with extension
Prius	CM500	✓	✓	✓	✓		
Prius	Elite	✓	✓	✓	✓		Foot extension*
Select Patient Care	Contour Stretcher	X	X	✓	X		
Steris	Horizon Series	X	X	✓	X		Column bed
Stralus	C 200 Series Hospital Acute Care Bed	✓	✓	✓	✓		
Stryker	Eye Surgery Stretcher	X	X	✓	X		Column bed
Stryker	S3	X	X	✓	X		Column bed
Unique Care	Standard High/Low Bed	✓	✓	✓	✓		
Volker	S962-2	✓	✓	✓	✓		



*Foot extension: Operators must lock foot extension before hitching, and ensure hitches are locked tight and secure from the widest point of the bed frame before operating.



Phone: 1300 934 471 (within Australia) | International: +61 3 9300 8500
 sales@fallshaw.com.au | www.fallshawgroup.com.au

MOBILITY SOLUTIONS TO HELP YOU
 PUSH, PULL & LIFT WITH EASE

