

Gzunda GZS PATIENT MOVER



The Gzunda GZS patient mover allows one person to effortlessly move heavy hospital beds, stretchers, treatment chairs or trolleys. The motorised unit operates equally well with or without a patient, simplifying the task of moving beds or patients between wards, departments, or to and from theatre, allowing the driver to focus on the needs of the patient during transport.

With a capacity of 600 kg, the silent motorised unit is able to move a hospital bed up or down ramps, over carpet, into lifts, and around tight corners with minimal effort by the operator. The GZS has easy-to-use controls, including key start, twist grip variable speed control, horn, and emergency stop button. As with all Electrodrive equipment, the GZ30 comes with a comprehensive 12-month warranty, including parts and labour.



Typical applications

Suitable for moving hospital beds, stretchers, treatment chairs or trolleys. No under-bed clearance is required.

Features

Load capacity	Move up to 600 kg
Max. speed	Up to 5 km/hour
Speed mode	Three speed control with forward, reverse and emergency stop.
Usability	<ul style="list-style-type: none"> • Highly manoeuvrable in tight spaces. • Intuitive controls with twist grip throttle and horn. • No under-bed clearance required. • No driver's licence required. • Simple to use. • Quiet, smooth operation. • Zero emissions. • Non-marking tyres.
Hitching	<ul style="list-style-type: none"> • Push-button hitching—silent, fast and efficient. • A variety of hitches are available.
Dimensions (L/W/H)	<ul style="list-style-type: none"> • GZSIN(NH): 1245/920/1230 mm • GZSUS(NH): 1245/1015/1230 mm
Battery	Two 12V 33Ah MK-gel batteries with 24V smart charger.

Safety features

- Emergency stop button.
- Emergency back-off button.

ORDER CODES

Gzunda GZS 600 kg split tine — no hitch (international)	GZSIN(NH)
Gzunda GZS 600 kg split tine — no hitch US)	GZSUS(NH)